



# Apolipoprotein A1 Antigen (ApoA1 Ag) (DAG-WT116)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunogen/Calibrators/Control
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements.
Size	1 mg
Buffer	20mM PBS, pH7.2
Preservative	0.05% Proclin300
Storage	2-8 °C (up to 6 Months), below -20 °C (up to 2 Years). Avoid repeat freeze thaw cycles.

## BACKGROUND

Introduction	Apolipoprotein A1 is a protein that in humans is encoded by the APOA1 gene. s the major component of HDL particles, it has a specific role in lipid metabolism. The text in a 2014 report suggested that APOA1 mRNA is regulated by endogenously expressed antisense RNA.
--------------	---