



Cefazolin [BSA] (DAG-WT1160)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Applications	LF, ELISA
Format	Liquid
Size	1 mg
Buffer	PBS buffer
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Cefazolin is a first-generation cephalosporin antibiotic used for the treatment of a number of bacterial infections. Specifically it is used to treat cellulitis, urinary tract infections, pneumonia, endocarditis, joint infection, and biliary tract infections. It is also used to prevent group B streptococcal disease around the time of delivery and before surgery. It is typically given by injection into a muscle or vein.
Keywords	Cefazolin; Cephalosporin antibiotic; Cefacidal