



Recombinant RSV prefusion glycoprotein F0 protein (DAG-WT1165)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90 % as determined by SDS-PAGE
Conjugate	His
Applications	ELISA, CLIA, CG
Molecular Weight	56.1 kDa
Format	Liquid
Size	1 mg
Buffer	PBS, pH 7.4
Preservative	None
Storage	Store at -80°C
Ship	Dry ice

BACKGROUND

Introduction	Human respiratory syncytial virus (RSV) is a virus that causes respiratory tract infections. It is a major cause of lower respiratory tract infections and hospital visits during infancy and childhood. A prophylactic medication (not a vaccine) exists for preterm (under 35 weeks gestation) infants, infants with certain congenital heart defects (CHD) or bronchopulmonary dysplasia (BPD), and infants with congenital malformations of the airway. Treatment is limited to supportive care (for example C-PAP), including oxygen therapy.
--------------	--

Keywords

Fusion protein; Respiratory syncytial virus; RSV; preF
