



Recombinant Krebs Von den Lungen-6 Antigen (KL-6 Ag) (DAG-WT117)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	His
Applications	Calibration, Control
Format	Liquid
Size	1 mg
Buffer	20mM PB, 150mM NaCl, pH 7.4
Preservative	None
Storage	2-8 °C (up to 1 Month), below -20 °C (up to 18 Months). Avoid repeat freeze thaw cycles.

BACKGROUND

Introduction	Krebs von den Lungen-6 (KL-6) is an indicator of lung inflammation or damage. Elevated levels are associated with various types of interstitial lung diseases (ILD).
Keywords	Krebs von den Lungen-6; KL-6