



Human Bence Jones Kappa protein (DAG-WT1171)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 95 % as determined by SDS-PAGE
Conjugate	N/A
Applications	ID, ELISA, WB, DIBA
Reconstitution	Reconstituted by adding 0.5 ml of sterile distilled water, mixed gently by inversion until complete dissolution is obtained.
Format	Lyophilized
Size	1 mg
Preservative	None
Storage	After arrival store at 4°C. Prolonged storage may be at -20°C. Avoid multiple freeze/thaw cycles.
Ship	Wet ice

BACKGROUND

Introduction	Bence Jones protein is a type of protein that is produced by plasma cells, which are a type of white blood cell. When present in high levels in urine, it can be a sign of multiple myeloma, a type of blood cancer. The term "human Bence Jones" refers to this protein when it is found in humans. Overall, the presence of human Bence Jones protein can be an important diagnostic tool for identifying multiple myeloma, and ongoing research continues to explore its potential use in developing new treatments and therapies for this disease.
---------------------	--

Keywords

Bence Jones protein; BJP; Myeloma
