



## Cefapirin sodium (DAG-WT1174)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Purity	> 98%
Conjugate	N/A
Cas.No	24356-60-3
Applications	Standard
Molecular Weight	445.445 g/mol
Format	Solid
Size	100 mg
Storage	Store at 2°C to 8°C
Ship	Wet ice

### BACKGROUND

Introduction	Cephalosporins are beta-lactam antimicrobials used to manage a wide range of infections from gram-positive and gram-negative bacteria. The five generations of cephalosporins are useful against skin infection, resistant bacteria, meningitis, and other infections.
Keywords	Cephalosporin; Beta-lactam antimicrobials