



Recombinant Mouse Alpha-synuclein protein (DAG-WT1175)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	Unconjugated
Applications	ELISA, LF
Molecular Weight	15.4 kDa
Format	Lyophilized
Size	100 µg, 1 mg
Buffer	Lyophilized from PBS
Preservative	None
Storage	Store at 2 to 8°C for one week. Store at -20 to -80°C for twelve months from the date of receipt.
Ship	Wet ice

BACKGROUND

Introduction	Alpha-synuclein is a protein that, in humans, is encoded by the SNCA gene. Alpha-synuclein is a neuronal protein that regulates synaptic vesicle trafficking and subsequent neurotransmitter release.
Keywords	Alpha-synuclein; SNCA; alpha SYN

GENE INFORMATION

UniProt ID	O55042
------------	--------
