



Recombinant of Human Neutrophil Gelatinase-Associated Lipocalin (NGAL) protein (DAG-WT118)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|----------------------|---|
| Purity | > 97%, as determined by SDS-PAGE |
| Conjugate | Unconjugated |
| Applications | Calibrator or standard in immunoassays |
| Format | Liquid |
| Concentration | Batch dependent - please inquire should you have specific requirements. |
| Size | 1 mg |
| Buffer | 0.01 M PBS, pH 7.4 |
| Preservative | None |
| Storage | Store at -80°C. Avoid freeze/thaw cycles |

BACKGROUND

| | |
|---------------------|---|
| Introduction | Neutrophil Gelatinase-Associated Lipocalin Neutrophil gelatinase-associated lipocalin (NGAL) is protein bound to gelatinase, and it was first described in neutrophils. |
|---------------------|---|