



Recombinant HDV sHDAg [His] (Cap) (DAG-WT1191)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA, CLIA
Format	Liquid
Size	1 mg
Buffer	20mM PBS
Preservative	None
Storage	Store at -20°C to -80°C

BACKGROUND

Introduction	Hepatitis D occurs only among people who are infected with the Hepatitis B virus. Transmission requires contact with infectious blood. At-risk populations include intravenous drug abusers and people who have received multiple blood transfusions. Symptoms include abdominal pain, nausea, and fatigue.
Keywords	Hepatitis D virus; HDV