



Recombinant canine thyroid-stimulating hormone (TSH) (DAG-WT1204)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Canine
Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Format	Liquid
Size	100 µg
Buffer	Tris/PBS-based buffer
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Thyroid-stimulating hormone (also known as TSH or thyrotropin) is a peptide hormone synthesized and secreted by thyrotrope cells in the anterior pituitary gland, which regulates the endocrine function of the thyroid gland.
Keywords	Thyroid stimulating hormone; TSH