



Recombinant Feline Troponin I (DAG-WT1208)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Feline
Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Format	Liquid
Size	1 mg
Buffer	Tris/PBS-based buffer
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Troponin I is a protein that is used as a diagnostic marker for myocardial infarction or heart attack. It is a subunit of the troponin complex, which is found in cardiac and skeletal muscle and is involved in muscle contraction. When heart muscle cells are damaged, troponin I is released into the bloodstream, where it can be measured by a blood test.
Keywords	Feline; Troponin I; cTnI

GENE INFORMATION

UniProt ID	Q863B6
------------	--------