



Recombinant Uricase (UA) [His] (34kDa) (DAG-WT121)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	His
Applications	Immunoassay
Molecular Weight	34kDa
Bio-activity	10.0 U/mg
Format	Lyophilized Powder
Size	1 mg
Storage	Store at -20°C. Avoid freeze/thaw cycles

BACKGROUND

Introduction	The enzyme urate oxidase (UO), uricase or factor-independent urate hydroxylase, absent in humans, catalyzes the oxidation of uric acid to 5-hydroxyisourate
--------------	---