



Recombinant PINP Antigen (DAG-WT1210)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	18 kDa
Format	Liquid
Size	100 µg, 500 µg, 1mg
Buffer	25mM Tris-HCl, 150mM NaCl, 50%Glycerol,pH 8.0
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Amino-terminal propeptide of type I procollagen (PINP) is released into blood circulation during bone formation. PINP is used as a bone turnover marker for the assessment of fracture risk and monitoring of osteoporosis treatment. PINP is recommended as reference bone formation marker by IOF and IFCC.
Keywords	Amino-terminal propeptide of type I procollagen; PINP