



# Native Krebs Von den Lungen-6 Antigen (KL-6) (DAG-WT1216)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|              |                |
|--------------|----------------|
| Conjugate    | N/A            |
| Applications | TIA            |
| Format       | Liquid         |
| Size         | 5KU, 10KU      |
| Buffer       | PBS            |
| Preservative | None           |
| Storage      | Store at -20°C |

## BACKGROUND

|              |                                                                                                                                                                      |
|--------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Introduction | Krebs von den Lungen-6 (KL-6) is an indicator of lung inflammation or damage. Elevated levels are associated with various types of interstitial lung diseases (ILD). |
| Keywords     | Krebs von den Lungen-6; KL-6                                                                                                                                         |