



Recombinant Human GPC3 VLP (DAG-WT1232)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by HPLC
Conjugate	N/A
Applications	ELISA, SPR
Molecular Weight	14 kDa
Format	Liquid
Size	100 µg, 500 µg
Buffer	20mM HEPES, 500mM NaCl, 10% Glycerol, 0.1% Tween20, pH 7.7
Preservative	None
Storage	Stable at -80°C for 12 months

BACKGROUND

Introduction	Glypican-3 (GPC3), a member of the glypican family that attaches to the cell surface by a glycosylphosphatidylinositol anchor, is overexpressed in HCC cases and is elevated in the serum of a large proportion of patients with HCC. GPC3 expression contributes to HCC growth and metastasis.
Keywords	GPC3; Glypican 3

GENE INFORMATION

