



Recombinant Human CD24 VLP (DAG-WT1239)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by HPLC
Conjugate	N/A
Applications	ELISA, SPR
Molecular Weight	5 kDa
Reconstitution	Reconstituting to a concentration more than 100 ug/ml is recommended. Dissolve the lyophilized protein in distilled water.
Format	Lyophilized
Size	100 µg, 500 µg
Buffer	Lyophilized from 0.22um filtered solution in PBS, pH 7.4. Normally 8% trehalose is added as protectant before lyophilization
Preservative	None
Storage	Reconstituted protein stable at -80°C for 12 months, 4°C for 1 week

BACKGROUND

Introduction	This gene encodes a sialoglycoprotein that is expressed on mature granulocytes and B cells and modulates growth and differentiation signals to these cells. The precursor protein is cleaved to a short 32 amino acid mature peptide which is anchored via a glycosyl phosphatidylinositol (GPI) link to the cell surface. This gene was missing from previous genome assemblies, but is properly located on chromosome 6. Non-transcribed pseudogenes have
--------------	---

been designated on chromosomes 1, 15, 20, and Y. Alternative splicing results in multiple transcript variants.

Keywords	CD24; CD24 Molecule
-----------------	---------------------

GENE INFORMATION

UniProt ID	P25063-1
-------------------	----------
