



Native Human Thyroid Peroxidase (DAG-WT1247)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 95%
Conjugate	N/A
Applications	Immunoassays
Reconstitution	Reconstitute in distilled water
Format	Liquid
Size	100 µg
Buffer	20mM Tris buffer, 50mM KCl, pH 7.5
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Thyroid peroxidase (TPO) is an enzyme normally found in the thyroid gland. TPO plays an important role in the production of thyroid hormones. A TPO test detects antibodies against TPO in the blood.
Keywords	Thyroid peroxidase; TPO