



## Native AFP L3 antigen (DAG-WT1249)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

Conjugate	N/A
Applications	Calibration, Standard
Format	Liquid
Size	10 µg, 50 µg, 100 µg
Buffer	Protein stabilizer
Preservative	None
Storage	Store at -20°C. Avoid multiple freeze/thaw cycles.

### BACKGROUND

Introduction	AFP-L3 is a subfraction of AFP that is produced by malignant hepatocytes. The ratio of AFP-L3 to total AFP is reported as a percentage, and over 10% AFP-L3 is consistent with a diagnosis of HCC.
Keywords	AFP-L3; Alpha-fetoprotein; AFP