



11-Dehydrothromboxane B2 [BSA] (DAG-WT1255)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	BSA
Applications	ELISA
Format	Liquid
Size	1 mg
Buffer	PBS buffer
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	11-Dehydrothromboxane B2 (or 11-dehydro-TXB2) is produced from the breakdown of thromboxane A2. It is released by activated platelets and urine levels of 11-dehydro-TXB2 can be used to monitor the response to aspirin therapy when used to prevent heart disease and in diseases where platelet activation is prominent.
Keywords	11-Dehydrothromboxane B2; 11-dehydro-TXB2