



Native Human Aspartate Aminotransferase (AST) (DAG-WT1259)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	N/A
Molecular Weight	92 kDa
Format	Powder
Size	100 U, 1000 U
Storage	Store at -80°C up to 1 year

BACKGROUND

Introduction	AST (aspartate aminotransferase) is an enzyme that is found mostly in the liver, but it's also in muscles and other organs in your body. When cells that contain AST are damaged, they release the AST into your blood. An AST blood test measures the amount of AST in your blood.
Keywords	Aspartate aminotransferase; AST

GENE INFORMATION

UniProt ID	P00505
------------	--------