



Native 1-5-Anhydroglucitol-6-Phosphate Dehydrogenase (DAG-WT1262)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	One unit is defined as the amount of enzyme which converts 1 μ mole of 1,5 –AG–6–phosphate to C6H11O8P per minute at 37°C under the conditions specified in the assay procedure
Conjugate	N/A
Applications	Enzymatic determination
Molecular Weight	40 kDa
Format	Lyophilized powder
Size	1 KU
Buffer	Enzyme dilution buffer
Storage	Store at -80°C

BACKGROUND

Keywords 1,5-Anhydroglucitol-6-phosphate; 1,5-AG