



Recombinant DNase B (DAG-WT1280)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	N/A
Applications	Enzymatic determination
Format	Liquid
Size	1 KU
Buffer	Enzyme dilution buffer
Storage	Store at -20°C

BACKGROUND

Introduction	Anti-DNase B (antideoxyribonuclease-B antibody) is one of the most common of several antibodies that are produced by the body's immune system in response to an infection with group A Streptococcus (strep infection). The anti-DNase B test measures the amount of this antibody in the blood. It is typically done with or following an antistreptolysin O (ASO) test, another test to detect antibody to a streptococcal antigen.
Keywords	DNase B; Deoxyribonuclease-B