



Recombinant Human TSHR [Fc] (DAG-WT1333)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	hFc
Applications	CLIA
Molecular Weight	100 kDa
Format	Liquid
Size	200 µg, 1 mg
Buffer	20mM PBS, pH7.4, 20% glycerol
Preservative	None
Storage	Store at -20°C

BACKGROUND

Keywords	TSHR; Thyroid stimulating hormone receptor
----------	--