



Purified bovine casein (DAG-WT1336)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	N/A
Applications	ELISA
Stability	≥ 12 months
Format	Liquid
Size	0.5mg, 1 mg
Buffer	Tris-glycine buffer
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Casein is a family of related phosphoproteins (α S1, aS2, β , κ) that are commonly found in mammalian milk, comprising about 80% of the proteins in cow's milk and between 20% and 60% of the proteins in human milk. Sheep and Cow milk have a higher casein content than other types of milk with human milk having a particularly low casein content.
Keywords	Casein; Dairy milk