



Recombinant ICA (Islet Cell Autoantigen) (DAG-WT1344)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA, CLIA, CG
Format	Liquid
Size	1 mg
Buffer	20 mM Tris-HCl, pH 8.0
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Islet cell autoantibodies and what they tell us. Islet autoantibodies are markers that appear when insulin producing beta cells in pancreas are damaged. They can be detected a long time before beta cells are completely destroyed. We use autoantibodies to estimate an individual's risk of developing type 1 diabetes.
Keywords	Islet cell autoantigen; ICA