



# Recombinant Human ADGRD2 VLP (DAG-WT1385)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|               |                                                                        |
|---------------|------------------------------------------------------------------------|
| Species       | Human                                                                  |
| Purity        | > 85% as determined by SDS-PAGE                                        |
| Conjugate     | N/A                                                                    |
| Applications  | N/A                                                                    |
| Stability     | 12 months                                                              |
| Endotoxin     | < 1 EU/μg                                                              |
| Format        | Liquid                                                                 |
| Concentration | Batch dependent - please inquire should you have specific requirements |
| Size          | 100 μg, 1 mg                                                           |
| Storage       | Store at -80°C                                                         |

## BACKGROUND

|          |                                           |
|----------|-------------------------------------------|
| Keywords | ADGRD2; VLP; Multi-pass transmembrane VLP |
|----------|-------------------------------------------|

## GENE INFORMATION

|            |        |
|------------|--------|
| UniProt ID | Q7Z7M1 |
|------------|--------|