



Recombinant Human ADGRG7 VLP (DAG-WT1401)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 85% as determined by SDS-PAGE
Conjugate	N/A
Applications	N/A
Stability	12 months
Endotoxin	< 1 EU/μg
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	100 μg, 1 mg
Storage	Store at -80°C

BACKGROUND

Keywords	ADGRG7; VLP; Multi-pass transmembrane VLP
----------	-------------------------------------------

GENE INFORMATION

UniProt ID	Q96K78
------------	--------