



Recombinant Human C5AR2 VLP (DAG-WT1428)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 85% as determined by SDS-PAGE
Conjugate	N/A
Applications	N/A
Stability	12 months
Endotoxin	< 1 EU/μg
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	100 μg, 1 mg
Storage	Store at -80°C

BACKGROUND

Keywords	C5AR2; VLP; Multi-pass transmembrane VLP
----------	--

GENE INFORMATION

UniProt ID	Q9P296
------------	--------