



Human FSHR full length protein-synthetic nanodisc (DAG-WT1493N)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Flag
Applications	ELISA
Molecular Weight	78.3 kDa
Format	Lyophilized
Size	50 µg, 100 µg
Buffer	Lyophilized from 20 mM Tris-HCl, 150 mM NaCl, pH 8.0
Storage	Store at -80°C

BACKGROUND

Keywords	FSHR; Follicle stimulating hormone receptor
----------	---

GENE INFORMATION

UniProt ID	P23945
------------	--------