



CTX-I Antigen (DAG-WT1767)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 95%
Conjugate	N/A
Applications	ELISA
Molecular Weight	1942.1 Da
Format	Liquid
Size	10 µg
Buffer	PBS
Storage	Store at -20°C

BACKGROUND

Introduction	CTX-I can be used as an indicator of bone resorption status as well as an aid in monitoring bone resorption changes during hormone replacement and bisphosphonate therapies.
Keywords	CTX-I; C-terminal telopeptide I