



Recombinant *Dermatophagoides farinae* Der f 2 (DAG-WT1793)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	≥ 95% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	15 kDa
Format	Lyophilized
Size	250 µg, 1 mg
Buffer	Lyophilized from sterile PBS, pH 7.4
Preservative	None
Storage	Store at -20°C to -80°C
Ship	Wet ice

BACKGROUND

Keywords	<i>Dermatophagoides farinae</i> ; Der f 2
----------	---

GENE INFORMATION

UniProt ID	Q00855.1
------------	----------