



Recombinant *Felis domesticus* Fel d 1 (DAG-WT1796)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	≥ 90% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	21.6 kDa
Format	Lyophilized
Size	250 µg, 1 mg
Buffer	Lyophilized from sterile 20mM Tris, 50mM NaCl, 10% glycerol, pH 7.4
Preservative	None
Storage	Store at -20°C to -80°C
Ship	Wet ice

BACKGROUND

Keywords	<i>Felis domesticus</i> ; Fel d 1
----------	-----------------------------------