



Recombinant *Betula pendula* Bet v 1 (DAG-WT1798)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	≥ 90% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	17.6 kDa
Format	Lyophilized
Size	250 µg, 1 mg
Buffer	Lyophilized from sterile 50 mM Tris, pH 8.0
Preservative	None
Storage	Store at -20°C to -80°C
Ship	Wet ice

BACKGROUND

Keywords	<i>Betula pendula</i> ; Bet v 1
----------	---------------------------------

GENE INFORMATION

UniProt ID	P15494.2
------------	----------