



# Recombinant *Alternaria alternata* Alt a 1 (DAG-WT1799)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA
Format	Lyophilized
Size	250 µg, 1 mg
Buffer	Lyophilized from sterile 20 mM Tris, 300 mM NaCl, pH 7.5
Storage	Store at -20°C to -80°C

## BACKGROUND

Keywords	<i>Alternaria alternata</i> ; Alt a 1
----------	---------------------------------------