



# Recombinant *Aspergillus fumigatus* Asp f 1 (DAG-WT1801)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Purity	≥ 90% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	18.3 kDa
Format	Lyophilized
Size	250 µg, 1 mg
Buffer	Lyophilized from sterile PBS, pH 7.4
Preservative	None
Storage	Store at -20°C to -80°C
Ship	Wet ice

## BACKGROUND

Keywords	<i>Aspergillus fumigatus</i> ; Asp f 1
----------	--

## GENE INFORMATION

UniProt ID	CAA41217.1
------------	------------