



Recombinant *Penaeus aztecus* Pen a 1 (DAG-WT1804)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	≥ 95% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	36 kDa
Format	Lyophilized
Size	250 µg, 1 mg
Buffer	Lyophilized from sterile 20mM Tris, 50mM NaCl, 2M Urea, pH7.4
Preservative	None
Storage	Store at -20°C to -80°C
Ship	Wet ice

BACKGROUND

Keywords	<i>Penaeus aztecus</i> ; Pen a 1
----------	----------------------------------