



Recombinant Triticum aestivum Tri a 14 (DAG-WT1805)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	≥ 95% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	9 kDa
Format	Lyophilized
Size	250 µg, 1 mg
Buffer	Lyophilized from sterile 20mM Tris, 500mM NaCl, 30% glycerol, 5mM β2-SH, pH 7.5
Preservative	None
Storage	Store at -20°C to -80°C
Ship	Wet ice

BACKGROUND

Keywords	Triticum aestivum; Tri a 14
----------	-----------------------------