



Recombinant *Blomia tropicalis* Blo t 2 (DAG-WT1806)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	≥ 95% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	15 kDa
Format	Lyophilized
Size	250 µg, 1 mg
Buffer	Lyophilized from sterile PBS, pH 7.4
Preservative	None
Storage	Store at -20°C to -80°C
Ship	Wet ice

BACKGROUND

Keywords	<i>Blomia tropicalis</i> ; Blo t 2
----------	------------------------------------

GENE INFORMATION

UniProt ID	AAQ73483.1
------------	------------