



Recombinant *Blomia tropicalis* Blo t 21 (DAG-WT1808)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	≥ 95% as determined by SDS-PAGE
--------	---------------------------------

Conjugate	His
-----------	-----

Applications	Immunoassays
--------------	--------------

Molecular Weight	16.3 kDa
------------------	----------

Format	Lyophilized
--------	-------------

Size	250 µg, 1 mg
------	--------------

Buffer	Lyophilized from sterile PBS, pH 7.4
--------	--------------------------------------

Preservative	None
--------------	------

Storage	Store at -20°C to -80°C
---------	-------------------------

Ship	Wet ice
------	---------

BACKGROUND

Keywords	<i>Blomia tropicalis</i> ; Blo t 21
----------	-------------------------------------

GENE INFORMATION

UniProt ID	AAX34047.1
------------	------------