



# Recombinant Dermatophagoides pteronyssinus Der p 5 (DAG-WT1809)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Purity	≥ 95% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	16.6 kDa
Format	Lyophilized
Size	250 µg, 1 mg
Buffer	Lyophilized from sterile PBS, pH 7.4
Preservative	None
Storage	Store at -20°C to -80°C
Ship	Wet ice

## BACKGROUND

Keywords	Dermatophagoides pteronyssinus; Der p 5
----------	---

## GENE INFORMATION

UniProt ID	XP_027199623.1
------------	----------------