



Recombinant *Artemisia vulgaris* Art v 4.02 (DAG-WT1812)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	≥ 95% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	15.3 kDa
Format	Lyophilized
Size	250 µg, 1 mg
Buffer	Lyophilized from sterile 20mM Tris, 100mM NaCl
Preservative	None
Storage	Store at -20°C to -80°C
Ship	Wet ice

BACKGROUND

Keywords	<i>Artemisia vulgaris</i> ; Art v 4.02
----------	--