



Recombinant Phleum pratense Phl p 5 (DAG-WT1815)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	≥ 95% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	30.9 kDa
Format	Lyophilized
Size	250 µg, 1 mg
Buffer	Lyophilized from sterile PBS, pH 7.4
Preservative	None
Storage	Store at -20°C to -80°C
Ship	Wet ice

BACKGROUND

Keywords	Phleum pratense; Phl p 5
----------	--------------------------

GENE INFORMATION

UniProt ID	Q40960-1
------------	----------