



Recombinant *Actinidia deliciosa* Act d 8 (DAG-WT1819)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	≥ 90% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	26.7 kDa
Format	Lyophilized
Size	250 µg, 1 mg
Buffer	Lyophilized from sterile 20mM Tris, 300mM NaCl, pH 8.0
Preservative	None
Storage	Store at -20°C to -80°C
Ship	Wet ice

BACKGROUND

Keywords	<i>Actinidia deliciosa</i> ; Act d 8
----------	--------------------------------------