



Native Angiotensinogen (DAG-WT182)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	1st 12 Amino acids: Asp-Arg-Val-Tyr-Ile-His-Pro-Phe-His-Leu-Val-Ile (453 total)
Purity	> 95% (SDS-PAGE)
Conjugate	Unconjugated
Applications	Immunoassays
Molecular Weight	61-65 kDa
Format	Lyophilized
Size	100 µg
Preservative	None
Storage	Store at -20°C for long term. Avoid freeze / thaw cycles.

BACKGROUND

Introduction	Angiotensin is a peptide hormone that causes vasoconstriction and an increase in blood pressure. It is part of the renin–angiotensin system, which regulates blood pressure. Angiotensin also stimulates the release of aldosterone from the adrenal cortex to promote sodium retention by the kidneys.
Keywords	AGT; ANHU; SERPINA8; hFLT1; Angiotensinogen

GENE INFORMATION

Gene Name	AGT, ANHU, SERPINA8
-----------	---------------------

