



# Recombinant *Arachis hypogaea* Ara h 6.0101 (DAG-WT1821)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Purity	≥ 95% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	16.2 kDa
Format	Lyophilized
Size	250 µg, 1 mg
Buffer	Lyophilized from sterile 0.1% TFA, 30% ACN
Preservative	None
Storage	Store at -20°C to -80°C
Ship	Wet ice

## BACKGROUND

Keywords	<i>Arachis hypogaea</i> ; Ara h 6.0101
----------	--

## GENE INFORMATION

UniProt ID	Q647G9.1
------------	----------