



Recombinant *Carya illinoensis* Car i 1 (DAG-WT1822)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	≥ 85% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	59.5 kDa
Format	Lyophilized
Size	250 µg, 1 mg
Buffer	Lyophilized from sterile 20mM Tris, 300mM NaCl, pH 7.4
Preservative	None
Storage	Store at -20°C to -80°C
Ship	Wet ice

BACKGROUND

Keywords	<i>Carya illinoensis</i> ; Car i 1
----------	------------------------------------