



Recombinant *Vespula vulgaris* Ves v 5 (DAG-WT1827)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	≥ 95% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	27.1 kDa
Format	Lyophilized
Size	250 µg, 1 mg
Buffer	Lyophilized from sterile 20mMTris, 300mMNaCl, 10% glycerol, pH 7.5
Preservative	None
Storage	Store at -20°C to -80°C
Ship	Wet ice

BACKGROUND

Keywords	<i>Vespula vulgaris</i> ; Ves v 5
----------	-----------------------------------