



Recombinant *Blattella germanica* Bla g 4 (DAG-WT1829)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA
Molecular Weight	30 kDa
Format	Lyophilized
Size	250 µg, 1 mg
Buffer	Lyophilized from sterile 20 mM Tris, 140 mM NaCl, pH 8.0
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords	<i>Blattella germanica</i> ; Bla g 4
----------	--------------------------------------