



# Native Antithrombin III (ATIII) (DAG-WT183)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Purity	> 95% (SDS-PAGE)
Conjugate	Unconjugated
Applications	Immunoassays
Molecular Weight	58,000
Format	Lyophilized
Size	1 mg, 10 mg, 100 mg
Buffer	Lyophilized powder from Tris-HCl, pH 8.0 with 150 mM sodium chloride
Preservative	None
Storage	Store at -20°C for long term. Avoid freeze / thaw cycles.

## BACKGROUND

Introduction	Antithrombin is also termed Antithrombin III (AT III). The designations Antithrombin I through to Antithrombin IV originate in early studies carried out in the 1950s by Seegers, Johnson and Fell.
Keywords	Antithrombin III; AT III; Antithrombin

## GENE INFORMATION

Gene Name	SERPINC1, AT3, PRO0309
-----------	------------------------