



Recombinant Human Serum Amyloid A [His] (DAG-WT1831)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	His
Applications	ELISA, LFIA
Format	Lyophilized
Size	1 mg
Buffer	Lyophilized from sterile 20mM Tris,500mM NaCl,pH8.0
Preservative	None
Storage	Store at 2 to 8°C
Ship	Wet ice

BACKGROUND

Introduction	Serum amyloid A (SAA) proteins are a family of apolipoproteins associated with high-density lipoprotein (HDL) in plasma. Different isoforms of SAA are expressed constitutively (constitutive SAAs) at different levels or in response to inflammatory stimuli (acute phase SAAs). These proteins are produced predominantly by the liver. The conservation of these proteins throughout invertebrates and vertebrates suggests that SAAs play a highly essential role in all animals.
Keywords	Serum amyloid A; SAA